


AVANTIK CSF²⁰

AUTOMATIC FILM COVERSGLIPPER

The Avantik CSF20 Automatic Film Coverslipper is designed to provide high quality, coverslipped slides in an easy-to-use, automated solution.



1. Xylene Pump System
2. Starter Basket Holder for Input Slides
3. Output Tray for Empty Slide Basket
4. Waste Bottle
5. Output Tray for Finished Slides Basket
6. End Basket Loading Station
7. Coverslipping Film Holder
8. User-Friendly LCD Display

Product Features

- ❖ Coverslips up to 1,200 slides per hour
- ❖ No glues or mounting media needed
- ❖ Ideal for H&E, IHC and cytology slide coverslipping
- ❖ Large, clear front provides protection against contamination
- ❖ Displays process status and remaining capacity for slides, xylene and coverfilm

User Benefits

- ❖ Integrated carbon filter for standalone use
- ❖ Intuitive menu navigation via touch display
- ❖ LED interior lighting for great interior visibility
- ❖ Rapid slide drying for immediate storage
- ❖ All components (such as coverslipping film) can be easily inserted and exchanged
- ❖ Minimal maintenance required, compared to traditional glass coverslippers
- ❖ Convenient windows in the doors allow for quick exchange of new slide baskets without interrupting coverslipping process

Peace of Mind,
Guaranteed Every Day.
Avantik, **A Higher Standard.**



Part Number	CSF20
Slide Basket Size	20 pcs
Coverslipping Time	~3 seconds per slide
Coverslipping Medium Capacity	Bottle of 500 ml
Coverslipping Film (width x length)	24 mm x 70 m
Coverslipping Length	40 - 60 mm
Waste Bottle Capacity	100 ml
Dimensions (width x depth x height)	730 x 310 x 545 mm
Weight	143 lbs
Power Supply	115 - 230 V / 50-60 Hz

Accessories	Quantity	Part Number
Tissue-Tek Film Coverslipping, 60 m	5/Case	SL6070
Pathotape Film Coverslipping, 70 m	1/Each	SL6400
Medite Film Coverslipping, 66 m	1/Each	SL6543
Cutting Blade	5/Box	SE3153
Charcoal Filter, Disposable	1/Each	SY6319
Basket for 20 Slides, Plastic	1/Each	GL5074
Waste Bottle, 100 ml	1/Each	GL5075
Xylene Bottle, 500 ml	1/Each	GL5076
Medite Stainer Compatible Hanger	1/Each	EQ6156
Epredia Gemini Stainer Compatible Hanger	1/Each	EQ6157
Leica Stainer Compatible Hanger	1/Each	EQ6158



SCAN NOW

to learn more about
the Avantik CSF20

Contact Us Today!

avantik-us.com

800.783.9424

