


AVANTIK TFB

TISSUE FLOATATION BATH



The **Avantik Tissue Floatation Bath** is an LED-illuminated, rectangular stretching bath for paraffin sections that has been tried and tested in the routine of everyday laboratory work.

Product Features



Heated Work Surface

Coated with a resistant plastic and offers enough space to store slides for dyeing.



Removeable Glass Bowl

Conveniently removable and dishwasher safe.



Temperature Regulation

Temperature sensor constantly measures exact temperature indicated by a digital display.

Contact Us Today!
avantik-us.com
800.783.9424

 **AVANTIK**
A HIGHER STANDARD


Product Specifications

Part Number	TFB23
Bath Capacity	max 1.5 liters
Temperature	+30°C bis +55°C
Dimensions Heated Front Area (width x depth)	370 x 75 mm
Dimensions Glass Tray (width x depth x height)	270 x 165 x 50 mm
Dimensions Overall (width x depth x height)	430 x 365 x 160 mm
Weight	ca. 7 kg
Power Supply	115 V / 50 - 60 Hz / 250 VA



SCAN NOW

to learn more about the
Avantik Tissue Floatation Bath



Contact Us Today!

avantik-us.com
800.783.9424

 **AVANTIK**
A HIGHER STANDARD