


AVANTIK AURA MICROTOME BLADES

THE IDEAL BLADE FOR ROUTINE SECTIONING



Avantik Aura microtome blades are ideal for routine sectioning and precise cuts when durability and consistency are important.

Product Details

- Three-faceted cutting edge for optimal sharpness
- Coated with a PTFE technology formula to improve the life of the blade over multiple cuts
- No pre-cleaning of blades required; ready to use upon dispensing
- Ideal for sectioning of both hard and soft tissue
- Made in the USA

Resilience Aura Low Profile (LP)

Item #: SE3145

Qty: 50 Blades/Pack

Resilience Aura High Profile (HP)

Item #: SE3144

Qty: 50 Blades/Pack

Peace of Mind.
Guaranteed Every Day
Avantik, **A Higher Standard.**



Avantik Aura Microtome Blades

The Ultimate Microtome Blade

Obtaining high-quality sections relies on high-quality microtome blades. Aura microtome blades offer increased consistency, impressive durability, and precise performance. As a result, these blades are ideal for routine sectioning when resilience and reliability are essential.



SCAN NOW

to learn more about Aura Blades
or request a sample pack.

Contact Us Today!
avantik-us.com
800.783.9424

MKTG-SS-P0002

Peace of Mind.
Guaranteed Every Day
Avantik, **A Higher Standard.**