

SLIDE & CASSETTE STORAGE



Avantik offers several stackable laboratory storage cabinets for various needs.

14 Drawer Metal Slide Storage Cabinet - Stackable

Stackable laboratory storage cabinets for microscope slides. Powder coated metal cabinets designed to be stackable. Base, lid and cabinet can be purchased separately. **All units designed to stack together.**

| Item | Description | Qty | Capacity |
|--------|--|------|-------------|
| ST7067 | 14 Drawer Slide Storage Cabinet 485mm x 480mm x 126mm | Each | 5600 Slides |
| ST7071 | 14 Drawer Cassette Storage Cabinet 485mm x 480mm x 91mm | Each | 900 Blocks |
| ST7069 | Base Unit - 485mm x 480mm x 76mm | Each | n/a |
| ST7068 | Top Cover - 485mm x 480mm x 26mm | Each | n/a |
| ST7072 | Spring Divider | Each | n/a |
| ST7073 | Magnetic Slide/Cassette Stop | Pk/5 | n/a |



6 Drawer Metal Slide Storage Cabinet & Base - Holds 4,500 Slides

Easy-slide drawers, capacity 4,500 standard slides. **Stackable with Paragon and Fisher cabinets.**

| Item | Color | Qty |
|--------|---------------|------|
| ST7078 | White Cabinet | Each |
| ST7079 | White Base | Each |

6 Drawer Metal Slide Storage Cabinet & Base - Holds 5,000 Slides

Six drawer cabinet holds up to 5000 standard 3" x 1" microscope slides. Drawers slide smoothly and quietly. Plastic drawer dividers can be removed to store objects up to 2" wide. Interlocking, stackable to 10 cabinets high. Includes sponge blocks spacers and identification labels. Base Model ST7077 has rubber feet to prevent scratching and slipping. 2 year warranty. **Stackable with the Boekel cabinets.**

| Item | Color | Qty |
|--------|---------------|------|
| ST7076 | Cabinet White | Each |
| ST7077 | Base White | Each |



Call us today.

800-783-9424
info@avantik-us.com
avantik-us.com

